

# **TK-Col TH (Truking Hard) IEX**

## **Ion Exchange Pre-Packed Column**

### **Product Manual**



# Content

1. Product Introduction .....	3
2. Technical parameters .....	4
3. Methods of use .....	6
4. Cleaning and regeneration .....	8
5. Sterilization and storage .....	9
6. Destruction and recycling .....	9
7. Ordering Information .....	10

## 1. Product Introduction

Ion exchange chromatography (IEC) is a highly effective method for separating and purifying biomolecules. This technique primarily relies on interactions between positive and negative charges, exploiting differences in the nature and magnitude of charges carried by various biomolecules under specific conditions to achieve separation. It features high loading capacity, excellent resolution, controllable conditions, and scalability, and has been widely applied in pharmaceuticals, chemical engineering, metallurgy, food processing, and other fields.

The TH IEX Series matrix is a highly cross-linked, extremely rigid agarose. Compared to conventional IEX TA FF, it exhibits higher dynamic loading capacity under high flow rate conditions, enabling adaptation to high-flow operations. Its physical and chemical properties remain stable throughout purification and cleaning/regeneration processes, with the medium's loading capacity and performance unaffected by external buffers. This makes it more suitable for large-scale industrial production.

- TH IEX possesses the following characteristics:
  - High flow rate, large volume processing capacity, short operating cycles, and improved yield
  - High-rigidity substrate with superior pressure resistance
- The TH IEX series primarily includes the following three types of media:
  - TH-Q Strong Anion Exchange Media
  - TH-SP Strong Cation Exchange Media
  - TH-S Strong Cation Exchange Media
  - TH-DEAE Weak Anion Exchange Media

TK-Col TH IEX pre-packed columns are ready-to-use ion exchange chromatography columns featuring TH IEX series media packed into TK-EC 1ml, TK-EC 4.9ml, TK-EC 5ml, and TK-EC 20ml empty chromatography columns. This eliminates the hassle of self-packing columns and the risk of suboptimal column efficiency. These pre-packed columns are widely used in laboratory process development and small-scale sample preparation, suitable for the separation and purification of biomolecules such as peptides, recombinant proteins, and antibodies. They feature the following characteristics:

- Ready to use
- Stable column bed volume
- High flow rate
- High dynamic binding capacity
- Excellent physicochemical tolerance

## 2. Technical parameters

**Table 1: Technical Parameters for Various Media**

medium	TH-Q	TH-DEAE	TH-S	TH-SP
Appearance	Spherical, white paste-like substance that separates into layers when left to stand.			
Base frame	Highly rigid agarose containing long-chain glucan			
Media Type	Strong anion	weak anion	Strong cation	Strong cation
functional group	Quaternary amine	Diethylaminoethyl	sulfonylmethyl	Sulfopropyl
Ion load	160~200 $\mu\text{mol Cl}^-/\text{ml}$	290~350 $\mu\text{mol Cl}^-/\text{ml}$	110~140 $\mu\text{mol H}^+/\text{ml}$	170~220 $\mu\text{mol H}^+/\text{mL}$
Average particle size <sup>+</sup>	90 $\mu\text{m}$			
Withstand pressure	0.5 MPa			
Operating pH Range	2-12	2-12	4-12	4-13
Chemical Stability	Common aqueous solutions: 1 M NaOH, 6 M hydrochloric acid guanidine, 30% isopropanol, 70% ethanol			
pH Stability	2~14 (CIP); 2~12 (working)	2~14 (CIP); 2~9 (working)	3-14 (CIP); 4-12 (working)	3~14 (CIP); 4~13 (working)
Temperature tolerance	Operating temperature: 4–40°C. Do not freeze. Autoclavable at 121°C (add 0.1 M salt) (Anion: NaCl; Cation: NaAc)			
payload	Different proteins vary significantly, with binding capacities ranging from a few milligrams to over a hundred milligrams per milliliter of medium—all within the normal range.			
Recommended flow rate	90~500cm/h			
Storage	2~30°C, 20% ethanol or 2% benzyl alcohol		2~30°C, 20% ethanol or 2% benzyl alcohol, 0.2M NaAc	

**Table 2: Technical Parameters of Pre-Packed Columns (Part Numbers See Last Page)**

Product Name	Resin	Prepacked column volume (ml)	Inner diameter × Column bed height mm × mm	Recommended flow rate ml/min	Withstand pressure	Storage	Screen aperture (μm)	
TK-Col Hard Q	TH-Q	1	7×25	<3	0.5MPa (5bar) (72.5psi)	20% ethanol	10	
		4.9	8×100	<4				
		5	16×25	<16				
		20	16×100	<16				
TK-Col Hard DEAE	TH-DEAE	1	7×25	<3				20% ethanol containing 0.2 M sodium acetate
		4.9	8×100	<4				
		5	16×25	<16				
		20	16×100	<16				
TK-Col Hard S	TH-S	1	7×25	<3		20% ethanol containing 0.2 M sodium acetate		
		4.9	8×100	<4				
		5	16×25	<16				
		20	16×100	<16				
TK-Col Hard SP	TH-SP	1	7×25	<3				20% ethanol containing 0.2 M sodium acetate
		4.9	8×100	<4				
		5	16×25	<16				
		20	16×100	<16				

### 3. Methods of use

- ◆ To prevent clogging of the chromatography column, all samples and buffers must be filtered through a 0.45μm membrane.
- ◆ To achieve optimal separation, avoid significant temperature differences between the buffer and the chromatography column.
- ◆ Store the chromatography column away from direct sunlight.
- ◆ Pre-packed columns may be used in chromatography refrigerators, but the flow rate should be appropriately reduced.

#### 3.1 Connect the chromatography column to the chromatography system

- Open the packaging and remove the chromatography column.
- Inspect the column for damage and check if it has dried out due to air ingress during shipping. If either occurs, contact your Chutian Microsphere sales representative immediately.
- Secure the column next to the chromatography system, ensuring the flow direction is correct.
- Start the chromatography system, ensuring all air bubbles are purged. Set the system alarm pressure to 0.3 MPa, then adjust and maintain a flow rate of 0.2 mL/min.
- After purging air bubbles from the system, open both end caps of the column and connect it while maintaining the low flow rate.

#### 3.2 Preparation of Chromatography Columns

- Rinse: Chromatography columns are stored in 20% ethanol or 2% benzyl alcohol during transport. First, flush out the storage solution with 2 column volumes of distilled water.
- Sterilization: For sample safety, it is recommended to rinse with 0.5M NaOH for 2 column volumes upon first use, followed by flushing with 2 column volumes of distilled water.

#### 3.3 Column equilibration

- Buffer: Select buffer salts whose buffer ions do not interact with the medium's ligands. For equilibration buffers, use low-salt (less than 5 mS/cm) and high/low pH buffers (typically: cation exchange media—1 pH unit below the target's isoelectric point; anion exchange media—1 pH unit above the target's isoelectric point) to facilitate binding while considering sample stability in the buffer. Elution buffers are typically prepared by adding high-concentration salts (e.g., 1M NaCl) to the equilibration buffer.
- Flush the column with the equilibration buffer at the recommended flow rate. Equilibration is complete when the pH and conductivity of the eluting buffer match those of the buffer before entering the column, typically requiring 2 – 5 column volumes.

**Table 3: Buffers Suitable for Anion Exchange Chromatography**

pH Scope	Buffer salt	concentration (mM)	Balanced Ions	pKa(25°C )
4.3-5.3	N-Methylpiperazine	20	Cl <sup>-</sup>	4.75
4.8-5.8	Piperazine	20	Cl <sup>-</sup> or HCOO <sup>-</sup>	5.33
5.5-6.5	L-Histidine	20	Cl <sup>-</sup>	6.04
6.0-7.0	bis-Tris	20	Cl <sup>-</sup>	6.48
6.2-7.2	bis-Tris propane	20	Cl <sup>-</sup>	6.65;

8.6-9.6				9.10
7.3-8.3	Triethanolamine	20	Cl <sup>-</sup> or CH <sub>3</sub> COO <sup>-</sup>	7.76
7.6-8.6	Tris	20	Cl <sup>-</sup>	8.07
8.0-9.0	N-Methyldiethanolamine	20	Cl <sup>-</sup>	8.52
8.0-9.0	N-Methyldiethanolamine	50	Cl <sup>-</sup> or CH <sub>3</sub> COO <sup>-</sup>	8.52
8.4-9.4	Diethanolamine	20(pH8.4)	Cl <sup>-</sup>	8.88
		50(pH8.8)		
8.4-9.4	Propane 1,3-Diamino	20	Cl <sup>-</sup>	8.88
9.0-10.0	Ethanolamine	20	Cl <sup>-</sup>	9.50
9.2-10.2	Piperazine	20	Cl <sup>-</sup>	9.73
10.0-11.0	Propane 1,3-Diamino	20	Cl <sup>-</sup>	10.55
10.6-11.6	Piperidine	20	Cl <sup>-</sup>	11.12

**Table 4: Buffers Suitable for Cation Exchange Chromatography**

pH Scope	Buffer salt	concentration (mM)	Balanced Ions	pKa(25°C )
1.4-2.4	Maleic acid	20	Na <sup>+</sup>	1.92
2.6-3.6	Methyl malonic acid	20	Na <sup>+</sup> or Li <sup>+</sup>	3.07
2.6-3.6	Citric acid	20	Na <sup>+</sup>	3.13
3.3-4.3	Lactic acid	50	Na <sup>+</sup>	3.86
3.3-4.3	Formic acid	50	Na <sup>+</sup> or Li <sup>+</sup>	3.75
3.7-4.7 5.1-6.1	Succinic acid	50	Na <sup>+</sup>	4.21; 5.64
4.3-5.3	Acetic acid	50	Na <sup>+</sup> or Li <sup>+</sup>	4.75
5.2-6.2	Methyl malonic acid	50	Na <sup>+</sup> or Li <sup>+</sup>	5.76
5.6-6.6	MES	50	Na <sup>+</sup> or Li <sup>+</sup>	6.27
6.7-7.7	Phosphate	50	Na <sup>+</sup>	7.20
7.0-8.0	HEPES	50	Na <sup>+</sup> or Li <sup>+</sup>	7.56
7.8-8.8	BICINE	50	Na <sup>+</sup>	8.33

### 3.4 Flow velocity

- Based on the type of chromatography column, select a flow rate within the recommended range, as shown in Table 2.

### 3.5 Sample submission

- To prevent sample clogging of the column, samples must be filtered through a 0.45μm microporous membrane prior to injection. Adjust the sample's pH and conductivity to match the equilibration buffer (this can be achieved through dilution, ultrafiltration, or buffer replacement using TD-G25). Determine the injection volume based on the sample's substance content and the ion exchange resin's binding capacity.

### 3.6 Rinse

- Rinse with balanced buffer until the UV absorbance value decreases to an appropriate level.

### 3.7 Elution

- Elution: Linear gradient or stepwise gradient elution can be employed to increase the elution strength in the eluent, thereby eluting substances with varying binding strengths from the chromatographic column. Different components are collected, and the position of the target analyte is detected.

### 3.8 Rebirth

- Rinse the chromatography column with a high-concentration salt solution (e.g., 2M NaCl).
- Rebalance: After rinsing with the equilibration buffer, proceed with the second sample loading. Repeat as necessary.

### 3.9 Column Efficiency Evaluation

Column efficiency determination may employ acetone or NaCl as indicators. Prepare the indicator solution and mobile phase according to the table below.

**Table 5: Column Efficiency Determination Method**

Method	Acetone Method for Column Efficiency Testing	NaCl Method for Column Efficiency Testing
Sample	1.0% (v/v) acetone aqueous solution	0.8M NaCl (dissolved in water)
Sample volume	1.0% column volume	1.0% column volume
Mobile phase	Water	0.4M NaCl aqueous solution
Flow velocity	30 cm/h	30 cm/h
Test Data	UV 280 nm	Electrical Conductivity

### 3.10 Calculate column efficiency

Calculate the theoretical plate height (HETP), theoretical number of plates (N), and asymmetry factor (As) based on the UV or conductivity curve using the following formula::

$$HETP=L/N$$

$$N=5.54(VR/Wh)^2$$

Where: VR= Residual Volume

Wh= Half-peak width

L= Column height

N= Theoretical plate count


The units for VR and Wh should be consistent.

$$As=b/a$$

Where:

a = First half-width at 10% peak height

b = Second half-width at 10% peak height



## 4. Cleaning and regeneration

As the number of uses of the chromatography medium increases, contaminants continue to accumulate on the column. Regular in-situ cleaning prevents this buildup and maintains stable operating conditions. Determine



 <b>楚天微球</b> TRUKING MICRO-SPHERE	<b>TruKing Micro-sphere Biotechnology Co.</b> <b>Product manual</b>	Instruction Manual No: 713 Edition number : 01 Effective date : 2025.01.01
--	--	--

the frequency of in-situ cleaning based on the degree of contamination of the medium (if contamination is severe, it is recommended to perform in-situ cleaning after each use to ensure reproducibility of results).

The following cleaning conditions are recommended for different types of impurities and contaminants:

- First flush out tightly bound proteins with 2 column volumes of purified water.
- For removal of strongly hydrophobic proteins and precipitated proteins: First wash with 2–3 column volumes of 1M NaOH, then immediately flush with 5–10 column volumes of purified water.
- For removal of lipoproteins and lipid substances: First wash with 5-10 column volumes of 70% ethanol or 30% isopropanol, followed by 5-10 column volumes of purified water.
- Alternatively, combine the above two washing conditions by washing with a 30% isopropanol solution containing 1M NaOH.

*Note: 70% ethanol or 30% isopropyl alcohol should undergo degassing prior to use; flow rates during in-situ cleaning may be selected between 30-60 cm/h; reverse flushing may be employed for severe blockages.*

## 5. Sterilization and storage

TK-Col TH IEX All pre-installed columns can be treated with 1M NaOH for 0.5–1 hour or longer to achieve sterilization and pyrogen removal.

TK-Col TH IEX SP/S Pre-packed columns are stored in 20% ethanol containing 0.2M sodium acetate (or 10mM NaOH). To prevent ethanol evaporation and microbial growth, used columns should be replaced with fresh 20% ethanol solution every three months.

Other TK-Col TH IEX pre-packed columns should be stored in 20% ethanol (or 10 mM NaOH). To prevent ethanol evaporation and microbial growth, used columns should be replaced with fresh 20% ethanol solution every three months.

## 6. Destruction and recycling

Due to the difficulty in naturally degrading the packing material within the TK-Col TH IEX series pre-packed columns, incineration is recommended to protect the environment.

## 7. Ordering Information

**Table 7: Part Numbers and Packaging**

Pseudolaric acid	product number	wrap
<b>TK-Col Hard Q</b>	Y6348	1×1ml
	Y6045	5×1ml
	Y6046	1×5ml
	Y6047	5×5ml
	Y604703	1×4.9ml
	Y604704	1×20ml
<b>TK-Col Hard DEAE</b>	Y6349	1×1ml
	Y6048	5×1ml
	Y6049	1×5ml
	Y6050	5×5ml
	Y605003	1×4.9ml
	Y605004	1×20ml
<b>TK-Col Hard S</b>	Y6350	1×1ml
	Y6051	5×1ml
	Y6052	1×5ml
	Y6053	5×5ml
	Y605303	1×4.9ml
	Y605304	1×20ml
<b>TK-Col Hard SP</b>	Y6541	1×1ml
	Y6542	5×1ml
	Y6543	1×5ml
	Y6544	5×5ml
	Y6545	1×4.9ml
	Y6546	1×20ml