



楚天微球
TRUKING MICRO-SPHERE

Truking Micro-sphere Biotechnology Co.
Product manual

Instruction Manual No.: 919
Edition number: 01
Effective date: 2025.01.01

TK-Col GST 4FF

TK-Col 16/10 GST 4FF

TK-Col 26/10 GST 4FF

Affinity Chromatography Prepacked Columns

Product Manual



楚天微球生物技术(长沙)有限公司
TRUKING MICRO-SPHERE BIOTECHNOLOGY (CHANGSHA) CO., LTD

Content

1. Product Introduction	3
2. Technical parameters	4
3. Methods of use	5
4. Cleaning and regeneration	7
5. Sterilization and storage	8
6. Destruction and recycling	8
7. Ordering Information	8

1. Product Introduction

TA-GST 4FF is a medium made by coupling glutathione to highly cross-linked agarose gel, which is specially used for the isolation and purification of Glutathione S-Transferase (GST) and GST fusion proteins, and GST tags are commonly used for the expression of fusion proteins in modern genetic engineering, which is conducive to soluble expression and activity maintenance of proteins. GST tags are commonly used in modern genetic engineering to express fusion proteins, which is beneficial to the soluble expression and activity maintenance of proteins. For glutathione S-transferase and its fusion proteins from different sources, the target proteins can be purified in one step by using this medium, which is characterized by high pressure resistance, high flow rate, and mild operating conditions, which is conducive to the preservation of protein activity.

TK-Col 16/10 GST 4FF pre-packed columns are ready-to-use affinity columns filled with TA-GST 4FF media in TK-EC 16/20 chromatography column; TK-Col 26/10 GST 4FF pre-packed columns are ready-to-use affinity columns filled with TA-GST 4FF media in TK-EC 26/20 chromatography column. This series of columns eliminates the hassle of loading the columns by the customer and the risk of poor column performance. This type of pre-packed columns is widely used for laboratory process development, small amount of sample preparation, and suitable for the separation and purification of biomolecules such as GST-tagged proteins. It has the following features:

- Ready-to-use
- Volume stabilization of the column bed
- Good physical and chemical resistance

2. Technical parameters

Table 1 TA-GST 4FF technical parameters

Appearance	White slurry, layered when placed
Base frame	4% agarose
Particle size distribution range	45~165μm
Functional groups	Glutathione with 10 atomic arms
Dynamic Binding Load	> 10mg GST/ml filler
Pressure resistance	0.3 MPa
Chemical stability	Stabilized in commonly used water-soluble buffers: 1M HAc, 70% ethanol, 6M guanidine hydrochloride (1 hour at room temperature)
pH stability	3~12
Storage	2~30°C, 20% ethanol or 2% benzyl alcohol
Pressure Flow Rate	~450cm/h (TK-EC16/10 H=5cm 25°C)
Recommended Flow Rate Range	Sampling flow rate: <100cm/h

Table 2: Technical parameters for each pre-assembled column (see end page for item number)

Product name	Prepacked resin	Prepacked column volume ml	Inner diameter×Column bed height mm×mm	Recommended flow rate + ml/min	Storage	Pressure resistance	Sieve plate aperture (μm)
TK-Col GST 4FF	TA-GST 4FF	1	7×25	<1	2-8°C, 20% ethanol or 2% benzyl alcohol (for international shipments).	0.3MPa (3bar)	10
		4.9	8×100	<3			
		5	16×25	<5			
		20	16×100	<5			
TK-Col 16/10 GST 4FF		19.1-21.1	16×100 (±5)	<5			
TK-Col 26/10 GST 4FF		50.4-55.7	26×100 (±5)	<13			

 楚天微球 TRUKING MICRO-SPHERE	Truking Micro-sphere Biotechnology Co. Product manual	Instruction Manual No.: 919 Edition number: 01 Effective date: 2025.01.01
--	--	---

3. Methods of use

- ◆ *TK-Col 16&26 series chromatography columns are made of glass and should be handled gently to prevent breaking or affecting the column efficiency.*
- ◆ *To avoid clogging the column, all samples and buffers need to be filtered through 0.45um membrane.*
- ◆ *In order to get a good separation effect, avoid too much temperature difference between the buffer and the column.*
- ◆ *Keep the column out of direct sunlight.*
- ◆ *Chromatography columns can be used in a chromatography cooler; but the flow rate needs to be reduced appropriately.*

3.1 Connecting the column to the chromatography system

- Open the package and take out the column
- Check whether the column is intact, and whether the column has been dried out during transportation, if any of the above situations occurs, please contact Chutian Microsphere sales representative in time.
- Fix the column next to the chromatography system and pay attention to the flow direction of the column.
- Start the chromatography system, make sure the air bubbles in the chromatography system are drained, and set the alarm pressure of chromatography system to 0.3MPa, then adjust and keep the flow rate running at 0.2ml/min.
- After the chromatography system is purged of air bubbles, open the upper and lower plugs of the chromatography column and connect the chromatography column under low flow rate operation.

3.2 Pretreatment of chromatography columns

- Rinse, the chromatography column is stored in 20% ethanol or 2% benzyl alcohol (for international transportation) during transportation, first rinse off the storage solution with 2 column volumes of distilled water.
- Sterilization, for sample safety, it is recommended to rinse 2 column volumes with 0.5M NaOH before rinsing 2 column volumes with distilled water for the first use.

3.3 Equilibration of Chromatographic Columns

- Use the recommended flow rate to flush the chromatographic column with equilibrium buffer, the pH and conductance of the buffer to be exported are consistent with that of the buffer before entering the column that is to say that the column is well equilibrated, and generally 2~5 column volumes are needed.
- buffer: generally use neutral buffer, such as 20 mM PB, 0.15 M NaCl, pH 7.3.

3.4 Flow rate

- According to the type of chromatographic column, generally choose the recommended flow rate within the range of flow rate, for the height of the column in the 10~15cm high case can be selected on the sample flow rate selection <100cm/h, equilibrium, cleaning, elution can be appropriate to increase the flow rate, but do not exceed the flow rate of 300cm/h. The larger the column height, the slower the flow

rate. (See Table 2)

3.5 Sampling

- **Sample and sample volume:** The pH and conductivity of the sample need to be adjusted to be consistent with the binding buffer, and in order to prevent the sample from clogging the column, the sample needs to be filtered with a 0.45 μm microporous filtration membrane before sampling, and the volume of the sample is determined according to the impurity content in the sample and the binding loading of the medium.

3.6 Rinse

- Rinse with equilibration buffer until the UV absorption value drops to the appropriate value.

3.7 Eluent

- Reduced glutathione is commonly used for elution, e.g., 50 mM Tris, 10 mM reduced glutathione, pH 8.0. 1-10 mM DTT can be added to the buffer, which can increase the purity of the target.

3.8 Regeneration and rebalancing

- **Regeneration:** Regular in-situ cleaning prevents the accumulation of contaminants and maintains a stable working condition. In-situ cleaning should be performed after each use to ensure reproducible results. The general regeneration method is as follows: two column volumes of high pH buffer (0.1 M Tris-HCl, 0.5 M NaCl, pH 8.5) and low pH buffer (0.1 M sodium acetate, 0.5 M NaCl, pH 4.5) are washed three times alternately; 10 column volumes of binding buffer equilibrate the chromatography column.
- **Re-equilibration:** rinsing with the equilibration buffer is sufficient for a second up-sampling, and so on.

3.9 Column Effectiveness Evaluation

Column efficiency can be determined by using acetone as indicator or NaCl as indicator, and the indicator solution and mobile phase are prepared according to the following table.

Table 4: Column efficiency determination methods

Methods	Acetone Method for Column Efficacy	Column Efficacy by NaCl Method
Sample	1.0% (v/v) acetone in water	0.8M NaCl (dissolved in water)
Sample volume	1.0% column volume	1.0% column volume
Mobile phase	Water	0.4M NaCl aqueous solution
Flow rate	30 cm/h	30 cm/h
Detection Data	UV 280 nm	Conductivity

3.10 Calculating Column Effect

Theoretical plate height (HETP), theoretical number of plates (N) and asymmetry factor (As) were calculated from the UV or conductivity curves with the following equations:

$$HETP = L/N$$

$$N = 5.54(V_R/W_h)^2$$

Where: V_R = retained volume

W_h = half peak width

L = column height


N = theoretical plate number

The units of V_R and W_h should be the same;

$$A_s = b/a$$

Where: a = first half peak width at 10% peak height

b = second half peak width at 10% peak height



3.11 Evaluation of results

$$h = HETP/d_{50v}$$

d_{50v} = median particle size volume distribution (cm)

The h -value calculated by the above formula is less than 3, and the asymmetry factor is 0.8~1.8 then it is judged to be qualified. For unsatisfactory column efficiency the reason needs to be analyzed and the column reloaded.

4. Cleaning and regeneration

As the number of times the chromatography medium is used increases, the accumulation of contaminants on the chromatography column also increases. Regular in-situ cleaning can effectively prevent the accumulation of contaminants and maintain the stable working condition of the chromatography medium. Customers can determine the frequency of in-situ cleaning according to the degree of contamination of the media during use (if the contamination is more serious, it is recommended to carry out in-situ cleaning after each use to ensure the reproducibility of the results).

Regeneration: 2 column volumes of high pH buffer (0.1M Tris-HCl, 0.5M NaCl, pH 8.5) and low pH buffer (0.1M sodium acetate, 0.5M NaCl, pH 4.5) were washed alternately three times; 10 column volumes of binding buffer to equilibrate the chromatography column.

If the medium is used for a period of time and the protein binding capacity decreases due to excessive impurities deposited on the surface, the medium needs to be cleaned as follows:

Precipitate or denatured material:

Wash with 2 column volumes of 6M guanidine hydrochloride, followed by 5 column volumes of equilibration buffer;

Hydrophobically bound substances:

Wash with 2 column volumes of 70% ethanol followed by 5 column volumes of equilibration buffer.

5. Sterilization and storage

Since 20% ethanol or 2% benzyl alcohol preservation solutions are not bacteriostatic or pyrogenic, it is recommended that TA-GST 4FF media can be treated with 70% ethanol for 12h before and during use to reduce the risk of microbial contamination.

TA-GST 4FF media are sold with 20% ethanol or 2% benzyl alcohol as preservation solution. After use, TA-GST 4FF should be stored in 20% ethanol in an airtight container at 2-30°C. To prevent ethanol evaporation and microbial growth, it is recommended that the preservation solution be replaced with fresh preservation solution every 3 months.

6. Destruction and recycling

Since TA-GST 4FF is difficult to degrade in nature, incineration is recommended for environmental protection.

7. Ordering Information

Table 5 Article number and packaging

Product	Item No.	Norm
TK-Col GST 4FF	Y6233	1×1ml
	Y6235	5×1ml
	Y6234	1×5ml
	Y6236	5×5ml
	Y623603	1×4.9ml
	Y623604	1×20ml
TK-Col 16/10 GST 4FF	Y6237	1pac.
TK-Col 26/10 GST 4FF	Y6238	1pac.